MEDICAL PRODUCTS

BAUERFEIND®



PRODUCT CATALOG
MEDICAL COMPRESSION STOCKINGS

2019

BAUERFEIND.CA



BAUERFEIND.
FOR A BETTER QUALITY OF LIFE.

MEDICAL BRAND **SINCE 1929**

CUSTOMER SERVICE

OFFICE LOCAL
(450) 629-9889

OFFICE TOLL FREE
1-877-629-9889

FAX LOCAL
(450) 629-9895

FAX TOLL FREE
1-877-629-9895

E-MAIL
info@galiengroup.com



CANADA
GALIEN
2017 Rue Cunard
H7S 2N1
Laval, Quebec



CONTENT

4 **Bauerfeind – by people for people**
6 **The three pillars of compression therapy**
7 **A systematic approach to compression**

8 **TREATMENT IN PHLEBOLOGICAL STANDARD CASES**
10 **Compression Sock TRAINING** and size chart
12 **Compression Sock PERFORMANCE** and size chart
14 **VenoTrain® discrétion** and size chart
16 **The new color system**
18 **VenoTrain® look**
20 **VenoTrain® micro**
22 Design and Batik Edition
24 Fashion Colors 2019
26 **VenoTrain® business**

28 **TREATMENT IN PHLEBOLOGICAL SPECIAL CASES**
30 **VenoTrain® soft**
32 Fashion Colors 2019
34 **VenoTrain® soft arm**
36 **VenoTrain® soft S** and size chart
38 **VenoTrain® impuls** and size chart
40 **VenoTrain® ulcertec** and size chart

42 **LYMPHOLOGICAL TREATMENT**

45 **ACCESSORIES**
46 **VenoTrain® glider**
47 **VenoTrain® glider plus**

BAUERFEIND.
BY PEOPLE FOR PEOPLE.
FOR A BETTER QUALITY OF LIFE.

Bauerfeind AG is one of the world's leading manufacturers of medical aids. With our products and treatment concepts, we enable people to be more active and thus make an important contribution toward maintaining and restoring health.


MEDICAL KNOW-HOW AND TECHNICAL EXPERTISE
Bauerfeind combines decades of expertise in technology and manufacturing with innovation: the tradition of a medium-sized, family-run business with a global mindset. Bauerfeind's specialists have worked very closely with physicians, scientists, and developers right from the start. The result: patented, award-winning products and services.

From the thread to the completed product, our compression stockings are developed and manufactured in Germany.

EFFECTIVE VARIETY
Our VenoTrain compression stockings cover all uses/indications in modern compression therapy, offering medically proven effectiveness and unique material properties.

Thanks to the wide range of colors available, the compression products blend seamlessly into their users' everyday lives and can be combined with various fashion styles. For prevention or the treatment of lymphedema: the indication-specific solutions include a wide variety of circular knit and flat knit designs that are available as custom-made or standard products.

QUALITY WITHOUT COMPROMISE
Bauerfeind products don't simply meet standards – they fulfill the most stringent requirements in terms of quality, effectiveness, and wearing comfort. Regular inspection procedures ensure that our products earn the seals and marks of successful testing.



MADE IN GERMANY


KESSLER
BAUERFEIND














THE THREE PILLARS OF COMPRESSION THERAPY
AN ALL-IN-ONE SOLUTION

Bauerfeind products meet particularly high standards, allowing you to provide your customers with the best possible treatment. VenoTrain compression stockings offer indication-specific solutions for all areas of treatment in modern compression therapy: whether for prevention, for phlebological standard or special cases, or for conservative treatment of lymphedema or lipedema. The circular knit or flat knit compression stockings offer medically proven effectiveness, unique material properties, and a patient-oriented modern design.



01
TREATMENT IN PHLEBOLOGICAL STANDARD CASES (page 8 – 27)

The phlebological standard case refers to **mild to moderate venous disorders**. It also includes situations in which there is a particular risk for the development of vein problems, such as during pregnancy or traveling. VenoTrain compression stockings improve venous activity and combat venous leg problems.



02
TREATMENT IN PHLEBOLOGICAL SPECIAL CASES (page 28 – 41)

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse. That is why the VenoTrain range offers compression stockings especially for patients with restricted mobility, dry skin, particularly weak veins, or open wounds.




03
LYMPHOLOGICAL TREATMENT* (page 42 – 43)

Lymphedema or lipedema patients are generally dependent on compression therapy for their entire lives. The effectiveness and wearing comfort of the stockings is therefore particularly important, significantly impacting the patient's quality of life.
VenoTrain **flat knit and circular knit products** provide the required therapeutic compression, while offering a unique wearing comfort, and help the patient live an active life.

* Detailed information about products for lymphedema and lipedema can be found in the separate Lymphology Product Catalog.

SYSTEMATIC APPROACH TO COMPRESSION



* Partsch, H.: Klassifikation von chronischen Venenerkrankungen heute, Vasomed 1 / 2003
** Rabe, E. et al.: Bonner Venenstudie der DGP, Phlebologie 1 / 2003



TREATMENT IN PHLEBOLOGICAL STANDARD CASES
FOR MILD TO MODERATE VEIN PROBLEMS

The phlebological standard case refers to **mild to moderate venous disorders**. It also includes situations in which there is a particular risk for developing vein problems, such as during pregnancy or traveling. VenoTrain compression stockings have been designed precisely to improve venous activity and combat venous leg problems.



Color	Product Code	Length	Size
Coal/Rivera	**291800530000**	**2** = Short	**1** = Small
Silver/Rivera	**291800531000**	**3** = Long	**2** = Medium
Coal/Polar	**291800532000**		**3** = Large
Silver/Polar	**291800533000**		**4** = Xlarge
Coal/Pink	**291800197000**		
Silver/Pink	**291800196000**		

COMPRESSION SOCK **TRAINING**

MUSCLE TUNING MADE EASY!
Bauerfeind's Compression Sock Training supports your calf muscle and lower leg, and stimulates the parts of the muscles that are under particular stress. In addition to providing special benefits, the socks help to prevent premature fatigue and muscle cramping.
The result: a lower risk of injury and more energy throughout your activities.

DETAILS

TOE closed
MATERIAL CONTENT
93 % polyamide, 7 % elastane
CARE INSTRUCTIONS


FEATURES

- **Set pressure profile** – consistent therapeutic gradient (Ccl 1) for improved circulation
- **Lycra® sports material** – gentle on the skin and breathable
- **Taping zone** – more stability for the ankle thanks to a patented design
- **Precision fit** – anatomically correct for left and right leg
- **Muscle tuning zone** – for reduced muscle vibrations and thus less fatigue
- **Protective zone** – for less pressure and rubbing on the Achilles tendon
- **Shock-absorbing cushioning** – for greater comfort and protection

AVAILABLE IN



Length ℓD
short up to 41 cm (16")
long over 41 cm (16")
Type Knee high



CIRCUMFERENCE		S	M	L	XL
cC Calf	cm	31–40	34–44	37–47	40–51
	inch	12–15	13–17	14.5–18.5	15.5–20
cB Ankle	cm	20–22	22–25	25–28	28–31
	inch	8–9	9–10	10–11	11–12

C0	C1	C2	C3	C4	C5	C6



	Ccl	Color	Toe	Circumference	Length	Size
292 800 = Knee high	**1**	**10** = White	**00** = Closed toe	**0** = Normal	**2** = Short	**1** = Small
		17 = Black			**3** = Long	**2** = Medium
		34 = Rivera				**3** = Large
		35 = Red				**4** = Xlarge
		36 = Green				
		37 = Yellow				
		97 = Pink				

COMPRESSION SOCK **PERFORMANCE**

RUNNING MADE EASY!
The Bauerfeind Compression Sock Performance exerts noticeable compression on the leg that complies with medical standards (Ccl 1) and effectively improves circulation. This ensures an optimal blood supply to the muscles and helps them perform for longer, particularly in endurance disciplines such as cycling or marathons.

The result: More energy reserves for long distances and quicker regeneration.

DETAILS

TOE closed
MATERIAL CONTENT
75 % polyamide, 25 % elastane
CARE INSTRUCTIONS


FEATURES

- **Set pressure profile** – consistent therapeutic gradient (Ccl 1) for improved circulation
- **Calf stabilization** – reduces typical muscle vibrations and prevents premature fatigue
- **Breathable microfiber** – provides optimum thermoregulation and is hard-wearing and comfortably light to wear
- **Elegant design** – understated for professionals and goes perfectly with any style of clothing
- **Non-skid foot design** – comfortable protection from pressure points and during increased strain

AVAILABLE IN



Length ℓD
short up to 41 cm (16")
long over 41 cm (16")
Type Knee high



CIRCUMFERENCE		S	M	L	XL
cC Calf	cm	31–40	34–44	37–47	40–51
	inch	12–15	13–17	14.5–18.5	15.5–20
cB Ankle	cm	20–22	22–25	25–28	28–31
	inch	8–9	9–10	10–11	11–12

C0 C1 C2 C3 C4 C5 C6



STANDARD SIZE CHART: VenoTrain discrétion

LEG CIRCUMFERENCE

Size		1	2	3	4
cG	cm inch	42 – 56 16½ – 22	49 – 63 19 – 25	56 – 70 22 – 27½	63 – 77 29 – 30½
cC	cm inch	28 – 36 11 – 14	32 – 40 12½ – 15½	35 – 43 13½ – 17	38 – 46 15 – 18
cB	cm inch	18 – 21 7 – 8	21 – 24 8 – 9½	24 – 27 9½ – 10	27 – 31 10½ – 12

LENGTH

Size		short	long
AG | AT	cm inch	< 70 < 27½	> 70 > 27½
AD	cm inch	< 38 < 15	> 38 > 15

	Ccl	Color	Toe	Circumference	Length	Size
216 8002 = Knee high	**1**	**02** = Clair	**02** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
216 6072 = Thigh high plain band w/dots		**04** = Java/Nightshade			**3** = Long	**2** = Medium
216 1002 = Pantyhose		**11** = Beige				**3** = Large
		17 = Black				**4** = Xlarge
		22 = Espresso				

VenoTrain® discrétion
THE TRANSPARENT STOCKING

The exceptionally fine material of these compression stockings gives them a stylish transparency and lightness.

Available in various styles and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure, and everyday life.

COLORS



FEATURES

- Highly elastic, making the stockings easy to put on
- Hugs the legs without constricting them
- Breathable fine knit for a long-lasting wearing comfort

DETAILS

COMPRESSION CLASS
- CcI 1

TOE closed

MATERIAL CONTENT
- 70 % polyamide, 30 % elastane

CARE INSTRUCTIONS



STANDARD

	AD	AG	AT
Clair	●	●	●
Beige	●	●	●
Espresso	●	●	●
Black	●	●	●
Java	●	●	●



THE NEW COLOR SYSTEM
SIMPLE. CLEAR. COORDINATED.

SELECTING THE RIGHT COLOR HAS NEVER BEEN EASIER

The new VenoTrain color system makes it even easier for you and your customers to select the right color: Now, the circular knit VenoTrain medical compression stockings for mild to advanced vein problems – VenoTrain look, VenoTrain micro, VenoTrain soft, VenoTrain business and VenoTrain impuls – as well as the flat knit product VenoTrain curaflow, are all available in the colors Caramel, Navy, Anthracite and Black.

The popular variants VenoTrain micro and VenoTrain soft are also offered in the light skin tone Cream, the darker skin tone Amaretto, and the warm brown tone Espresso. If you love a fresh breeze, have a look at the limited edition for 2019 in ocean and coral.

	STANDARD COLORS				ADDITIONAL COLORS			FASHION COLORS	
VenoTrain®	Caramel	Navy	Anthracite	Black	Cream	Amaretto	Espresso	Ocean	Coral
									
look	●	●	●	●					
micro	●	●	●	●	●	●	●	●	●
business	●	●	●	●					
soft	●	●	●	●	●	●	●	●	●
impuls	●	●	●	●					
curaflow	●	●	●	●	●				

THE BENEFITS OF THE NEW COLOR SYSTEM FOR YOU:

- Coordinated color selection for the entire VenoTrain range
- Easy and quick overview of all products and their corresponding colors
- Personalized and straightforward customer advice



	Ccl	Color	Toe	Circumference	Length	Size
217 800 = Knee high	**1** = 20–30	**12** = Caramel	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
217 606 = Thigh high lace band	**2** = 30–40	**15** = Navy		**3** = Plus	**3** = Long	**2** = Medium
217 100 = Pantyhose		**17** = Black				**3** = Large
		69 = Anthracite				**4** = Xlarge

VenoTrain® look
ATTRACTIVE AND SHIMMERING

The VenoTrain look is a real highlight among medical compression stockings thanks to its glossy look – for a stylish appearance on any occasion!

The VenoTrain look enhances the figure with a push-up effect. Stretch comfort and its special trilobal fiber guarantee a perfect fit.

STANDARD COLORS

CARAMEL
NAVY
ANTHRACITE
BLACK

FEATURES

- Fashionable and elegant shimmer effect
- Slightly transparent

DETAILS

COMPRESSION CLASSES
- CcI 1
- CcI 2

TOE closed


MATERIAL CONTENT
- CcI 1: 60 % polyamide, 40 % elastane
- CcI 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS

STANDARD

AD AG AT

SIZE CHART



CIRCUMFERENCE		S NORMAL	S PLUS	M NORMAL	M PLUS	L NORMAL	L PLUS	XL NORMAL	XL PLUS
cG Thigh	cm	40–50	50–60	50–60	60–68	60–68	68–74	68–74	74–82
	inch	15.5–19.5	19.5–23.5	19.5–23.5	23.5–26.5	23.5–26.5	26.5–29.5	26.5–29.5	29.5–32.5
cC Calf	cm	28–33	33–38	33–38	38–43	38–43	43–48	43–48	48–54
	inch	11–13	13–15	13–15	15–17	15–17	17–19	17–19	19–21
cB Ankle	cm	19–22		22–25		25–28		28–32	
	inch	7.5–8.5		8.5–10		10–11		11–12.5	

C0 C1 C2 C3 C4 C5 C6



	Ccl	Color	Toe	Circumference	Length	Size
218 800 = Knee high	**1**	**11** = Cream	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
218 807 = Knee high w/silicone band 5 cm*	**2**	**12** = Caramel	**20** = Open toe	**3** = Plus	**3** = Long	**2** = Medium
218 605 = Thigh high lace band w/dots		**14** = Amaretto				**3** = Large
218 607 = Thigh high plain band w/dots**		**15** = Navy				**4** = Xlarge
218 100 = Pantyhose		**17** = Black				
218 130 = Maternity		**22** = Espresso				
		69 = Anthracite				
		Fashion Colors				

* Available only in black and caramel.
** Not available in fashion colors.

Tango	Ccl	Color	Toe	Circumference	Length	Size
284 800 = Knee high	**1**	Available in all colors	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
284 605 = Thigh high lace band w/dots	**2**			**3** = Plus	**3** = Long	**2** = Medium
284 100 = Pantyhose						**3** = Large
						4 = Xlarge

Jive	Ccl	Color	Toe	Circumference	Length	Size
285 800 = Knee high	**1**	Available in all colors	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
285 605 = Thigh high lace band w/dots	**2**			**3** = Plus	**3** = Long	**2** = Medium
285 100 = Pantyhose						**3** = Large
						4 = Xlarge

VenoTrain® micro
FASHIONABLE, MODERN, AND VERSATILE

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of our most popular compression stockings.

Available in many designs and a wide selection of colors, it highlights the individual style of your customers.

STANDARD COLORS

CARAMEL
ANTHRACITE
NAVY
BLACK

ADDITIONAL COLORS

CREAM
AMARETTO
ESPRESSO

FEATURES

- More than 50 per cent microfiber content: breathable, skin-friendly, and easy-care
- Figure-shaping panty area (pantyhose) and comfort sole
- Opaque full-fashioned stocking quality

DETAILS

COMPRESSION CLASSES
- CcI 1
- CcI 2

TOE closed, open


MATERIAL CONTENT
- CcI 1: 65 % polyamide, 35 % elastane
- CcI 2: 55 % polyamide, 45 % elastane

CARE INSTRUCTIONS

STANDARD

AD AG AT ATU

SIZE CHART



CIRCUMFERENCE		S NORMAL	S PLUS	M NORMAL	M PLUS	L NORMAL	L PLUS	XL NORMAL	XL PLUS
cG Thigh	cm	40–50	50–60	50–60	60–68	60–68	68–74	68–74	74–82
	inch	15.5–19.5	19.5–23.5	19.5–23.5	23.5–26.5	23.5–26.5	26.5–29.5	26.5–29.5	29.5–32.5
cC Calf	cm	28–33	33–38	33–38	38–43	38–43	43–48	43–48	48–54
	inch	11–13	13–15	13–15	15–17	15–17	17–19	17–19	19–21
cB Ankle	cm	19–22		22–25		25–28		28–32	
	inch	7.5–8.5		8.5–10		10–11		11–12.5	

C0 C1 C2 C3 C4 C5 C6


BATIK
NAVY




VenoTrain® micro
DESIGN EDITION AND BATIK

THE SECOND SKIN FOR TRENDSETTERS


A true fashion statement:
the Design Edition and individually dyed VenoTrain micro Batik versions create real highlights.

With a subtle Tango or Jive pattern, or wildly textured in batik with navy or black, your customers' legs will be real eye-catchers.






STANDARD COLORS




ADDITIONAL COLORS



STANDARD








VenoTrain® micro
FASHION COLORS 2019

COLOR WORKS WONDERS,
FOR YOU AND YOUR LEGS.

With fashion colors that are updated every two years, the sophisticated VenoTrain micro compression stockings make your legs look great and add that special something to your outfit. So it's easy to be fashion-conscious while taking care of your health.



COLORS

	OCEAN	CORAL
STANDARD		
TANGO		
JIVE		
BATIK NAVY		
BATIK BLACK		

STANDARD

AD AG AT







		Ccl	Color	Toe	Circumference	Length	Size
Long	**208 800** = Knee high	**1**	**12** = Caramel	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
Short	**208 810** = Knee high	**2**	**15** = Navy		**3** = Plus	**3** = Long	**2** = Medium
			69 = Anthracite				**3** = Large
			17 = Black				**4** = Xlarge

VenoTrain® business
TIMELESSLY ELEGANT

While traveling or at work: the VenoTrain business meets every requirement for a premium knee-high stocking. The combination of cotton and Timbrelle® high-tech fiber offers top-class wearing comfort with thermal regulation. Highly elastic LYCRA® and the extra-wide supportive band ensure that the stocking is easy to put on and stays perfectly in place.

STANDARD COLORS

CARAMEL
NAVY
ANTHRACITE
BLACK

FEATURES

- Innovative mix of materials
- Reinforced sole
- Timbrelle® high-tech fiber providing long-lasting freshness and healthy foot conditions

TIMBRELLE® CARE

DETAILS

COMPRESSION CLASSES
- CcI 1
- CcI 2

TOE closed


MATERIAL CONTENT
- CcI 1: 80 % polyamide, 14 % elastane, 6 % cotton
- CcI 2: 71 % polyamide, 24 % elastane, 5 % cotton

CARE INSTRUCTIONS

STANDARD

AD

SIZE CHART



CIRCUMFERENCE		S NORMAL	S PLUS	M NORMAL	M PLUS	L NORMAL	L PLUS	XL NORMAL	XL PLUS
cC Calf	cm	28–33	33–38	33–38	38–43	38–43	43–48	43–48	48–54
	inch	11–13	13–15	13–15	15–17	15–17	17–19	17–19	19–21
cB Ankle	cm	19–22		22–25		25–28		28–32	
	inch	7.5–8.5		8.5–10		10–11		11–12.5	

Available in shoe sizes 36–41 (short) or 41–49 (long)

C0 C1 C2 C3 C4 C5 C6



TREATMENT IN PHLEBOLOGICAL SPECIAL CASES
FOR MODERATE TO SEVERE VEIN PROBLEMS

Some patients with vein problems have **special requirements** that call for compression stockings with unique modes of action and properties. These requirements can be very diverse. That is why the VenoTrain special range offers compression stockings especially for patients with:

- **GENTLE AND EFFECTIVE – VenoTrain soft**
- **RESTRICTED MOBILITY – VenoTrain soft S**
- **PARTICULARLY WEAK VEINS – VenoTrain impuls**
- **OPEN WOUNDS – VenoTrain ulcertec**



	Ccl	Color	Toe	Circumference	Length	Size
212 800 = Knee high	**1**	**11** = Cream	**00** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
212 607 = Thigh high plain band w/dots	**2**	**12** = Caramel	**20** = Open toe	**3** = Plus	**3** = Long	**2** = Medium
212 620 = Thigh high w/waist attachment R		**15** = Navy				**3** = Large
212 630 = Thigh high w/waist attachment L		**17** = Black				**4** = Xlarge
212 100 = Pantyhose		**22** = Espresso				
		42= Anthracite				

VenoTrain® soft
GENTLE AND EFFECTIVE

The VenoTrain soft, always the right choice: this versatile stocking can withstand plenty of wear, yet it is very soft and gentle on the skin.

It is available in many styles with a proven medical effect and a beneficial massage action.

STANDARD COLORS

CARAMEL
ANTHRACITE
NAVY
BLACK

ADDITIONAL COLORS

CREAM
AMARETTO
ESPRESSO

FEATURES

- For men and women
- Gentle massage action
- Particularly gentle on the skin

DETAILS

COMPRESSION CLASSES
- CcI 1
- CcI 2

TOE closed, open


MATERIAL CONTENT
- CcI 1: 70 % polyamide, 30 % elastane
- CcI 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS

STANDARD

AD AG AGHB AT

SIZE CHART



CIRCUMFERENCE		S NORMAL	S PLUS	M NORMAL	M PLUS	L NORMAL	L PLUS	XL NORMAL	XL PLUS
cG Thigh	cm	40–50	50–60	50–60	60–68	60–68	68–74	68–74	74–82
	inch	15.5–19.5	19.5–23.5	19.5–23.5	23.5–26.5	23.5–26.5	26.5–29.5	26.5–29.5	29.5–32.5
cC Calf	cm	28–33	33–38	33–38	38–43	38–43	43–48	43–48	48–54
	inch	11–13	13–15	13–15	15–17	15–17	17–19	17–19	19–21
cB Ankle	cm	19–22		22–25		25–28		28–32	
	inch	7.5–8.5		8.5–10		10–11		11–12.5	

C0 C1 C2 C3 C4 C5 C6





VenoTrain® soft
FASHION COLORS 2019

SOFT
BUT WITH A
STRONG IMPACT.

If you love color, you won't be able to resist the fashion colours of the VenoTrain soft opaque compression stocking. The current range of shades are an elegant way to demonstrate that compression therapy can not only feel good, but look good too.

COLORS

OCEAN

CORAL

STANDARD

AD AG AGHB AT



Soft arm sleeve	Ccl	Color	Length	Size
212 554 = silicone band	**1**	**00** = Natural	**2** = Short	**1** = Small
	2	**12** = Caramel	**3** = Long	**2** = Medium
		17 = Black		**3** = Large
				4 = Xlarge

VenoTrain® soft Arm
EVERYDAY MOBILITY

Seamless with an extra-soft inner elbow area. This facilitates everyday mobility with complete medical efficacy.

COLORS



FEATURES

- Comfort elbow for greater freedom of movement, wearing comfort, and
- a reduced pressure sensation
- Noticeable effect throughout the entire wearing period
- Gentle on the skin, soft, and breathable

DETAILS

MATERIAL CONTENT
- Ccl 1: 74 % polyamide, 26 % elastane
- Ccl 2: 69 % polyamide, 31 % elastane

COMPRESSION CLASSES
- Ccl 1
- Ccl 2

CARE INSTRUCTIONS



STANDARD



CG

SIZE CHART

Length ℓG

short	39.0 – 44.0 cm (15 – 17")
long	44.1 – 50.0 cm (17 – 20")



		S min. – max.	M min. – max.	L min. – max.	XL min. – max.
cG	cm	24 – 30	27 – 33	29 – 36	32 – 40
	inch	9 – 12	11 – 13	11 – 14	13 – 16
cE	cm	22 – 26	24 – 29	27 – 32	30 – 35
	inch	9 – 10	9 – 11	11 – 13	12 – 14
cC	cm	14 – 15	15 – 16.5	16.5 – 18	18 – 20
	inch	5.5 – 6	6 – 6.5	6.5 – 7	7 – 8





	Ccl	Color	Toe	Circumference	Length	Size
212 800 = Knee high	**1**	**12** = Caramel	**04** = Closed toe	**2** = Normal	**2** = Short	**1** = Small
212 804 = Knee high w/silicone band (3 cm)	**2**		**24** = Open toe	**3** = Plus	**3** = Long	**2** = Medium
212 607 = Thigh high plain band w/dots						**3** = Large
						4 = Xlarge

VenoTrain® soft S
COMFORT IN THE EVENT OF RESTRICTED MOBILITY

Based on proven VenoTrain soft quality which offers a gentle massage effect to stimulate the blood flow, the VenoTrain soft S also takes into account special requirements of people who spend most of their time seated.
The wide heel and comfortable instep area ensure that the foot, which is angled during sitting, has enough space, and that no creases can develop.
The robust knitted fabric can withstand plenty of wear, and can also be easily put on by care staff.

STANDARD COLORS

Standard products only available in caramel

CARAMEL

BLACK

NAVY

ADDITIONAL COLORS

CREAM

ESPRESSO

FEATURES

- Wide heel area
- Adapted to the typical anatomical foot posture while sitting (90-degree angle)
- Breathable material with a revitalizing effect

DETAILS

COMPRESSION CLASSES
- Ccl 1
- Ccl 2

TOE closed, open

MATERIAL CONTENT
- Ccl 1: 70 % polyamide, 30 % elastane
- Ccl 2: 65 % polyamide, 35 % elastane

CARE INSTRUCTIONS


STANDARD

AD AG

SIZE CHART

Length ℓG
short up to 71 cm (28")
long over 71 cm (28")
Type Thigh high/pantyhose

Length ℓD
short up to 41 cm (16")
long over 41 cm (16")
Type Knee high



CIRCUMFERENCE		S NORMAL	S PLUS	M NORMAL	M PLUS	L NORMAL	L PLUS	XL NORMAL	XL PLUS
cG Thigh	cm	40–50	50–60	50–60	60–68	60–68	68–74	68–74	74–82
	inch	15.5–19.5	19.5–23.5	19.5–23.5	23.5–26.5	23.5–26.5	26.5–29.5	26.5–29.5	29.5–32.5
cC Calf	cm	28–33	33–38	33–38	38–43	38–43	43–48	43–48	48–54
	inch	11–13	13–15	13–15	15–17	15–17	17–19	17–19	19–21
cB Ankle	cm	19–22		22–25		25–28		28–32	
	inch	7.5–8.5		8.5–10		10–11		11–12.5	

C0 C1 C2 C3 C4 C5 C6




STANDARD SIZE CHART

Length ℓG

short	up to 70 cm (27.5")
long	over 70 cm (27.5")
Type	AG

Length ℓD

short	up to 39 cm (15.5")
long	over 39 cm (15.5")
Type	AD



CIRCUM-FERENCE		1		2		3		4		5		6	
		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
cG	cm	45–53	51–59	48–56	54–62	51–59	57–66	54–63	60–69	57–66	63–73	60–70	66–76
	inch	18–21	20–23	19–22	21–24	20–23	22–26	21–25	24–27	22–26	25–29	24–28	26–30
cC	cm	29–33	31–35	31–35	33–38	33–38	35–40	35–40	37–42	37–42	39–44	39–44	41–47
	inch	11–13	12–14	12–14	13–15	13–15	14–16	14–16	15–16.5	15–16.5	15–17	15–17	16–18.5
cB	cm	19–20.5		20.5–22		22–24		24–26		26–28.5		28.5–31	
	inch	7.5–8		8–9		9–9		9–10		10–11		11–12	

	Ccl	Color	Toe	Circumference	Length	Size
207 800 = Knee high	**2**	**12** = Caramel	**00** = Closed toe	**2** = Normal	**2** = Short	**1**
207 804 = Knee high w/silicone band (3 cm)	**3**	**17** = Black	**20** = Open toe	**3** = Plus	**3** = Long	**2**
207 607 = Thigh high plain band w/dots						**3**
207 620 = Thigh high w/waist attachment R						**4**
207 630 = Thigh high w/waist attachment L						**5**
						6

VenoTrain® impuls
POWERFUL SUPPORT FOR PARTICULARLY WEAK VEINS

Produces a particularly deep action by exerting high working pressure with even the slightest muscle movements. Easy to put on and breathable thanks to microfibers.

Highly effective for advanced vein problems, post-thrombotic syndrome, follow-up care for healed ulcers, and therapy in the early stages of lymphedema.

STANDARD COLORS

CARAMEL



BLACK



FEATURES

- The short-stretch knitted fabric generates a high working pressure with even the slightest muscle movements
- Skin-friendly microfiber material makes the product breathable and ensures temperature control

DETAILS

COMPRESSION CLASSES
- Ccl 2
- Ccl 3

TOE closed, open

MATERIAL CONTENT
- Ccl 2: 80 % polyamide, 20 % elastane
- Ccl 3: 75 % polyamide, 25 % elastane

CARE INSTRUCTIONS



STANDARD



AD **AG** **AGHB**


C0 C1 C2 C3 C4 C5 C6



STANDARD SIZE CHART

Length ℓD

short	up to 40 cm (16")
long	over 40 cm (16")
Type	AD



CIRCUM-FERENCE		XS		S		M		L		XL	
		NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS	NORMAL	PLUS
cC	cm	24–29	28–34	27–32	31–38	29–36	35–42	33–40	38–47	36–44	43–52
	inch	9–11	11–13	11–13	12–15	11–14	14–17	13–16	15–18.5	14–17.3	17–20
cB	cm	20–21.5		21.5–23		23–25		25–28		28–32	
	inch	8–8.5		8.5–9		9–10		10–11		11–13	

	Ccl	Color	Toe	Circumference	Length	Size
209 800 = Knee high w/liners	**0** = 30–40	**00** = Natural	**21** = Open toe (w/closed toe liners)	**2** = Normal	**2** = Short	**0** = Xsmall
				3 = Plus	**3** = Long	**1** = Small
						2 = Medium
209 800 = Liners	**0**	**10**	**01**			**3** = Large
						4 = Xlarge

VenoTrain® ulcertec
FOR VENOUS LEG ULCERS

The combination of an overstocking with a patented rhomboid knitted fabric and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.

The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

COLORS

NATURAL
OVERSTOCKING



WHITE
LINER

FEATURES

- Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort
- Differentiated treatment based on two designs within one compression class
- Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking


DETAILS

COMPRESSION CLASSES
- CcI 2

TOE
- Open (overstocking)
- Closed (liner)

MATERIAL CONTENT
- 87 % polyamide, 13 % elastane

CARE INSTRUCTIONS
- Overstocking:
- Liner:



STANDARD



AD

C0 C1 C2 C3 C4 C5 C6



LYMPHOLOGICAL TREATMENT

Bauerfeind has a wide range of compression products for lymphological conditions which take into account the patients' different requirements.



The circular knit **VenoTrain impuls** compression stocking combines a wide range of benefits: its breathable microfiber and firm compression thread offer unique wearing comfort with a high degree of effectiveness. Only very little movement is needed for the stocking to exert the required therapeutic pressure.



Thanks to our custom-made **VenoTrain curaflow** flat knit products, we can offer lymphatic patients perfectly tailored solutions with more than 1,000 different options of combining products and additional features. The beneficial therapeutic features of the range include the knitted fabric's high microfiber content, the flat sensitive stitching, and special comfort zones such as the extra-soft area behind the knee.

Comprehensive information and product details of our lymphology portfolio can be found in the Lymphology Product Catalog.

ACCESSORIES
THE PERFECT ADDITION FOR EVERYDAY USE

These handy extras from the VenoTrain accessories range are useful for ensuring that your patients feel comfortable using and wearing their VenoTrain compression stockings. They provide effective treatment support, make it easier to put on the products and care for the compression stockings.

01 Special detergent
cleans gently, ensuring careful treatment of your products, and providing lasting protection

02 Rubber gloves
help with putting the compression stockings on correctly



01



02

VenoTrain® glider
COMPACT AID FOR PUTTING ON AND TAKING OFF STOCKINGS

The VenoTrain glider is a practical aid for putting on and taking off medical compression stockings. It can be prescribed by a physician. The lightweight, hard-wearing material makes the stocking glide effortlessly over the arm or leg. It is suitable for compression stockings with an open or closed toe. Saving space and time, it is the ideal companion when you're out and about.

FEATURES

- Extremely lightweight, robust, and tear-resistant
- Easy to handle
- Ideal for travel

A video with donning instructions can be found at:
video.bauerfeind.com/venotrain-glider



PUTTING ON THE STOCKINGS: STEP 01



PUTTING ON THE STOCKINGS: STEP 02



PUTTING ON THE STOCKINGS: STEP 03



PUTTING ON THE STOCKINGS: STEP 04



PUTTING ON THE STOCKINGS: STEP 05



TAKING OFF THE STOCKINGS

VenoTrain® glider plus
ROBUST DONNING AID

The VenoTrain glider plus is a robust donning aid with an ergonomic design for all types of medical compression stockings. It can be prescribed by the physician.

Folding handles make it easy to put on the stockings, avoiding the need to bend down. It has been tested by an independent institute for ease of use, even for customers with restricted mobility.

FEATURES

- A compact, ergonomic design, light metal frame with folding handles, easy to use without the need for any great effort

Vous pouvez visionner une vidéo avec instructions d'enfilage sur :
video.bauerfeind.com/venotrain-glider-plus


PUTTING ON THE STOCKINGS: STEP 01


PUTTING ON THE STOCKINGS: STEP 02


PUTTING ON THE STOCKINGS: STEP 03


PUTTING ON THE STOCKINGS: STEP 04


PUTTING ON THE STOCKINGS: STEP 05


PUTTING ON THE STOCKINGS: STEP 06

BETTER QUALITY OF LIFE
FOR YOUR CUSTOMERS!

Have you got any questions, suggestions, or requests?
Please contact us:

CANADA
GALIEN
2017 Rue Cunard
H7S 2N1
Laval, Quebec

Office local (450) 629-9889
Office toll free +1 877 629-9889
Fax local +1 (450) 629-9895
Fax toll free +1 877 629-9895
info@galiengroup.com

Rev. 1 – 2019-05

